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# Content

<b>I. Introduction</b>	<b>22</b>
<b>II. eVito Body Weight Scale</b>	<b>23</b>
Display & functions	23
General information	24
Security information	26
<b>III. Setup &amp; Functions</b>	<b>27</b>
Changing the batteries	27
Setting up of the user data	28
Measuring	31
Wireless data transfer	32
Data transfer via user profile	32
Interpretation of measured values	33
<b>IV. Miscellaneous</b>	<b>37</b>
Maintenance	37
Warranty	38

# I. Introduction

Thank you that you have decided for the **eVito Body Weight Scale**. This product enables you to measure your body weight as well as content of adipose, content of water, content of muscle and bone mass. Please, read the instructions carefully before you first use the scale and keep the instructions for a future use.

## Principle of measurement

**eVito Body Weight Scale** is based on the principle of the bioelectric impedance analysis (BIA). A weak and non-perceptible and, therefore, not dangerous electric signal will be conducted through the whole body. The scale identifies the electric resistance.



Individual data such as age, body height, gender or degree of activity allows identifying the adipose content and the other parameters. Muscle tissue and water have a good electric conductivity and, therefore, a lower resistance than bone and adipose tissue, because adipose cells and bones hardly conduct current.

Please note, that the values of **eVito Body Weight Scale** are only an approximation to the real medical values of the body. Only a specialist who takes medical methods into account (e. g. computer tomography) can identify the exact parameters like content of adipose, content of water, content of muscles and bone mass.

**eVito Body Weight Scale** is an ANT+ device that supports the ANT+ device profile and it is compatible with receiving devices that support the ANT+. You can then transmit your values wirelessly to these receiving devices.

## II. eVito Body Weight Scale

### Display & functions

## Specifications

1. Device: **eVito Body Weight Scale**
2. Maximum capacity: 150 kg / 330 lb
3. Units: kg / lb (in 100 g / 0,2 lb steps)
4. Buttons: Set, Go up, Go down
5. Activity indices: 4 choices
6. Content of adipose: in % for 0,1 % steps
7. Content of water: in % for 0,1 % steps
8. Content of muscle mass: in % for 0,1 % steps
9. Basal metabolic rate: in kcal with 1 kcal per step
10. Active metabolic rate: in kcal with 1 kcal per step
11. User profiles: 8
12. Battery: 4 x AAA batteries (1,5 V)
13. Size (L x W x H): 325 x 325 x 19 mm
14. Weight: 1950 g
15. Operating temperature: 10°C - 40°C
16. Storage temperature: -20°C - 60°C
17. Symbol when overburdened: „----“
18. Symbol for low batteries: „Lo“
19. Communication/Data transfer: ANT+ wireless

## General Information

- Please place the scale on an even, hard floor. The floor is the basis for a correct measurement.
- If possible, please use the scale at the same time every day (we recommend in the morning), after having used the toilet, with empty stomach and without clothes in order to be able to compare the results.
- Please wait about 15 minutes after getting up before you measure your body weight so that the water of your body can disperse.
- The adipose can only be detected when you are barefoot. You can moisten the soles of your feet. Completely dry feet can result an incorrect measurement because they do not have enough conductivity.
- Please do not move during the measurement.
- Please wait several hours before measuring after unaccustomed physical exercise.

- The following persons or groups of persons will not get differing or not plausible results concerning their adipose and the other values:
  - children under 10 years
  - competitive athletes & bodybuilders
  - pregnant woman
  - people having fever, dialysis patients, people with edema symptoms or osteoporosis
  - persons taking cardiovascular medicine
  - persons taking vasodilatory or vasoconstrictive medicine
  - persons with considerable anatomic differences of their legs that can influence the height of the body (different lengths of the legs).
- The device is for your own use only and not for medical or commercial use.
- **eVito Body Weight Scale** corresponds to the EC directive 89/336.

## Security information

The **eVito Body Weight Scale** must not be used by persons with medical implants (e. g. cardiac pacemaker). Their function could be affected.

Do not use during pregnancy. The amniotic liquor inside the body can influence the measurements.

Please do not use the device with wet feet – slip hazard!

Do not swallow the batteries – danger of life. Please store batteries and device at place which children cannot reach. In case a battery has been swallowed, please consult a physician immediately.

Keep children away from the packaging materials - danger of asphyxiation!

Do not recharge the batteries. Do not take them to pieces. Do not throw them into a fire. Do not short-circuit them.

## III. Setup & Functions

### Changing the batteries

- Please open the cover of the battery case at the back of the device.
- Insert 4 AAA batteries (1,5 V).
- Please pay attention to the polarity.
- Close the cover of the battery case.

**Note:** In case the scale does not operate please remove the batteries and insert them again.

Your device has a symbol indicating that the batteries need to be changed. In case the batteries are almost empty, the symbol “Lo” can be seen on the display. You should change the batteries then.

Please remove the batteries in case you are not going to use the scale for one month or more in order to avoid an automatic discharge of the batteries.

The device as well as the batteries must be disposed according to the local regulations.

In order to guarantee a correct measurement you should calibrate the scale always after any location change. Just tip on the scale and wait until it shows “0.0” on the scale.

## Setting up of the user data

In order to measure the adipose content and other parameters, you have to insert your individual user data first. You can store up to 8 user profiles and recall the measured data. To set up the scale, please press “SET” (on the side of your scale) to create a new user profile.

### *User profile:*

The display now shows several illuminated preadjustments. The symbol “P1” appears on the display (bottom left) for user profile 1 (see Pict. 1). You can confirm this choice by pressing “SET” or you can create a new user profile by increasing the number (arrow keys at the side of the scale).

Pict. 1

### *Gender:*

Press „SET“ to confirm your choice and then chose your gender. The display shows the symbol for a man. If you want to change the gender please use one of the arrow keys until you can see the symbol for a woman. Press “SET” to confirm your choice.

### *Body height:*


Please insert your body height by using the arrow keys. You can insert body heights from 100 cm up to 250 cm. When you continue to press the buttons the indicated values go faster. Please confirm your value by pressing “SET”.

### *Age:*

The display shows the value for your age. Please use again the arrow keys to insert your age and confirm by pressing “SET”.

## Activity index:

Then the symbol for the choice of the activity index appears in the left display segment (see Pict. 2).

Pict. 2

Value defines a certain level of activity that is important for a medium term or a long term consideration. You can choose between 4 activity levels for each profile by using the arrow keys. Please confirm your choice and press "SET".

Defining the different levels of activity:

- Level 1: No physical exercise
- Level 2: Only little exercise such as walks, working in the garden or gymnastics.
- Level 3: Medium physical exercise (30 minutes physical exercise from 2 up to 4 times per week minimum)
- Level 4: High physical exercise (30 minutes physical exercise from 4 up to 6 times per week minimum)

Please note that the fat calculation is depending on the activity index. The higher the index the lower the fat calculation will be.

**ATTENTION:** When using your **eVito Body Weight Scale** with the **eVito system**, always choose Level 1 in order to have your evaluations in the **eVito system** working properly.

## Function setup

After the input of the activity index you will enter the function setup. The function setup can be set to ON or OFF. Switch through the arrows to ON or OFF and then confirm with „SET“.

**Heart Beat :** This is the automatic setting function.


Pict. 3

**Beep 1 :** A long beep will signal when the measurement is finished.


Pict. 4

**Beep 2 :** 3 short beeps for transferring only weight or 5 short beeps for transfer of all data. This indicates that the transmission was successful.

Pict. 5

**Beep 3 :** 2 long beeps will signal when the transmission was NOT successful. Switch through the arrows to ON or OFF and then confirm with „SET“.

Pict. 6

After having finished the setting up of your data the symbol “0,0” appears on the display and your **eVito Body Weight Scale** is ready to start the measurement.

## Measuring

After having determined all parameters and user profiles body weight, adipose and the other values can be measured.

- Please press “SET” and choose your user profile by using the arrow keys. Do not confirm by pressing “SET” since you would have to insert and confirm all the parameters again.
- The measurement starts automatically after a few seconds. First of all, **eVito Body Weight Scale** detects your body weight (left display segment). After that the measurement of the other parameters will start. You will see three circles on the display segments concerned. The scale is detecting the other parameters – content of water, content of muscles and bone mass (see Pict. 7).

- Finally all values are displayed (see Pict. 8).
- After a few seconds the indication of the body weight changes into the indication of the calorie consumption - basal metabolic rate (see Pict. 9), then into the indication of the active metabolic rate (in kcal). In addition the activity index is indicated by the symbol “activity index” (see Pict. 10).

## **Wireless data transfer**

Your **eVito Body Weight Scale** has an integrated ANT module which enables you to transfer measurement data to your computer, mobile phone or sport watch. There you can save, analyse and process the values in the **eVito system**.

### *What is ANT?*

ANT is a wireless data transfer standard. In order to use it, both transmission and receiving device must be equipped with an ANT module.

### **Data transfer via user profile**

Set up your user profile according to the instructions in this manual. The device is preset to transfer the data automatically. Please choose your user profile and start the measurement. Your device will transfer the measured data to the receiving device automatically after the measurement.

## Interpreting the results

The following guiding values help you to interpret the measured values of your **eVito Body Weight Scale**. These values are only a general recommendation; you should eventually discuss and interpret the data with your physician.

*Adipose content in %:* 

The following values are guiding values concerning your adipose content:

Age	Woman				Man			
	<i>under-nourished</i>	<i>healthy / normal</i>	<i>over-weigh</i>	<i>too much adipose</i>	<i>under nourished</i>	<i>healthy / normal</i>	<i>over-weigh</i>	<i>too much adipose</i>
10 - 12	< 12	12 - 22	22 - 30	> 30	< 8	8 - 18	18 - 24	> 24
13 - 18	< 15	15 - 25	25 - 33	> 33	< 8	8 - 18	18 - 24	> 24
19 - 30	< 20	20 - 29	29 - 36	> 36	< 8	8 - 18	18 - 24	> 24
31 - 40	< 22	22 - 31	31 - 38	> 38	< 11	11 - 20	20 - 26	> 26
41 - 50	< 24	24 - 33	33 - 40	> 40	< 13	13 - 22	22 - 28	> 28
51 - 60	< 26	26 - 35	35 - 42	> 42	< 15	15 - 24	24 - 30	> 30
60+	< 28	28 - 37	37 - 47	> 47	< 17	17 - 25	26 - 34	> 34

Table 1

## Content of water in your body in %:

The content of water inside your body is usually within the following ranges (only the normal range is indicated):

Age	Woman	Man
10 - 12	> 60	> 64
13 - 18	> 58.5	> 63.5
19 - 30	> 56	> 62.5
31 - 40	> 53	> 61
41 - 50	> 52	> 60
51 - 60	> 51	> 59
60+	> 50	> 58

Table 2

**Note:** Adipose has relatively low water content. For this reason the content of water in case of a high adipose % can be below the guiding values above. For the endurance athletes, the values can be higher than the guiding values because of the low adipose and the high content of muscles.

## Content of muscles in %:

The content of the muscles usually stays within the following ranges (only the normal range is indicated):

Age	Woman	Man
10 - 99	> 34	> 40

Table 3

## Bone mass in kg:

Our bones as well as our whole body are subject to the ageing process and all its consequences. During the childhood the bone mass is increasing fast and reaches its maximum at the age of 30-40. Then the bone mass is decreasing. You can influence the ageing process by a healthy nutrition (esp. calcium and vitamin D) and by a continuous physical exercise.

 Bone mass and bone density are different factors. Please do not confuse them. The latter can only be detected by a medical examination (e. g. computer tomography, ultra sound). For this reason the **eVito Body Weight Scale** cannot be used to identify changes of the bones or their durability. The bone mass can hardly be influenced but is changing according to factors like weight, body height, age or gender.

The bone mass can hardly be influenced but is changing according to factors like weight, body height, age or gender.

	Woman			Man		
Body weight	< 45 kg	45 - 60 kg	> 60 kg	< 60 kg	60 - 75 kg	> 75 kg
Bone mass	1,8 kg	2,2 kg	2,5 kg	2,5 kg	2,9 kg	3,2 kg

Table 4

### ***Basal Metabolic Rate (BMR):***

BMR means the amount of energy the body needs for its basal vital functions when the body is in the state of absolute calmness (e. g. when you stay in bed for 24 hours). This value mainly depends on factors like weight, body size and age.

This value is displayed on your **eVito Body Weight Scale**. The unit is kcal/day. The basis for the calculation is the science-based Harris-Benedict Formula. If your body gets less energy for a longer time, it can affect negatively your health.

### ***Active Metabolic Rate (AMR):***

The active metabolic rate means the amount of energy the body needs daily in the state of activity. The energy consumption increases with increased physical exercise. Your **eVito Body Weight Scale** has 4 different activity indices to determine the appropriate degree of activity.

In order to maintain the current weight the consumed energy must be replaced in form of food and drink. In case when this energy will not be replaced, the body uses the stored fat – the body loses weight. In case when the body gets more energy than the required AMR, the body is not able to use this energy – it will be stored inside the body and you gain weight.

### ***Time as a factor:***

Please note that only a long term result counts. Short term changes of your weight within a few days are absolutely normal and are often due to changes of the content of water. Only middle term or long term changes are related to the content of your adipose or your muscles. In case your weight decreases for a short time and the content of the adipose of your body increases or stays the same you only lost water (e.g. after a physical exercise or after having visited a sauna). If your weight increases in the medium term and your content of adipose decreases or stays the same you might have built up more muscle mass.

## **IV. Miscellaneous**

### **Maintenance**

- Please clean your device regularly by a damp cloth. Do not use caustic cleaning solutions.
- Do not put the device under running water.
- Protect the device from impacts, humidity, dust, chemicals or fluctuations of temperature. Do not store near heat sources such as ovens or heating installations.
- Never take the device to pieces by yourself. In this case the warranty will not be valid anymore. Only the manufacturer or the authorized distributor is allowed to repair the **eVito Body Weight Scale**. Please check the batteries before a complaint and change them if necessary.
- The **eVito Body Weight Scale** is in accordance with the EC directive 89/336. Please feel free to contact the manufacturer or the customer service for further questions regarding the applications of the scale.

The product may not be disposed of with the regular household garbage. Please dispose of the device according to the regulation for the waste of electrical and electronic equipment 2002/96/EC – WEEE.

## Warranty

HMM Diagnostics GmbH puts emphasis on a high quality standard. For this reason HMM Diagnostics GmbH gives a warranty of 2 years. Please store the warranty card included in your product set.

*Parts subject to regular wear and tear, batteries etc. are not subject to the warranty.*

## Manufacturer:

HMM Diagnostics GmbH  
Friedrichstr. 89  
D-69221 Dossenheim, Germany

[www.hmm.info](http://www.hmm.info)